



WHEELCHAIR CUSHIONS



Miracle Memory Cushion

Manufactured in USA by:

FLA ORTHOPEDICS®
Miramar, FL 33025


www.Posturainfo.com

Part # B5-125550

©2003 FLA Orthopedics, Inc., Miramar, FL 33025

Miracle Memory Cushion C1255 SERIES

*For Positioning,
Skin Protection,
& Pressure Relief*



CARE INSTRUCTIONS:

Spot wipe clean removable cover with damp cloth. Machine wash outer cover ONLY in cold water with mild soap on gentle cycle. Air dry away from heat.

FEATURES:

- **Slow Recovery Viscoelastic Foam** Enhances Comfort and Provides Superior Support - 45% More Foam by Density Than Other Foam Cushions
- **Revolutionary Micro-Cell Foam** Conforms to the Body's Contours and Absorbs Shock
- **Contoured Sitting Surface** for Proper Positioning, Increased Stability and Lower Extremity Control
- **Convex Seat Bottom** Contours a Sling Surface to Reduce "Hammocking"
- **Maximum Pressure Relief** Over High-Risk Areas
- **Low Profile**
- **Lightweight** - Easy to Handle and Transport
- **Fluid Proof Cover with Non-Slip Bottom** - Removable and Washable
- **Cover Meets CAL 117 Flammability Standard**
- **Carrying Handle**

INDICATIONS:

For extra seated comfort and support. For individuals at low or medium risk of developing pressure ulcers or who need or prefer slight leg positioning for proper posture and added comfort. For those individuals who use a wheelchair for mobility but can ambulate and shift weight independently.

CUSHION DIRECTION:

Place the cushion in the desired seat with the Medial Thigh Support to the front of the chair.

SIZE: Model Width x Depth x Thickness
C1255 18" x 16" x 2"

Cover Material: Top – 100% Treated Polyester with Polyurethane Fluid Proof Film
Bottom – Coated Polyester